

# FagronLab™ FG

Semi-automatic manual capsule filling system



## FagronLab™ FG

- A new cost-effective solution for quick and easy production of 100 capsules. Standard equipment includes a range of useful accessories such as powder tray, tamper tool, spreading spatulas
- Easy and quick capsule size change. Smartly expandable with purchase of extra low cost changing parts for the production of additional capsule sizes
- Flexibility; suitable for capsule sizes 00, 0, 1, 2, 3, 4
- Reduced powder static and minimised powder loss, high quality construction made of stainless steel
- Food grade certified plastic parts, compliant to EU Regulation 10/2011
- Over 99% capsule separation
- User friendly
- Warranty 2 years

**Output per hour**  
1.500 - 2.500 pcs

**Dimensions (HxLxW)**  
10x25x18cm

**Weight:** 2kg