


# FagronLab™ LF-500 UV

Positive pressure workbench



Positive Pressure Airflow Pattern

## FagronLab™ LF-500 UV

- HEPA Filter 99.999% efficiency at 0.3um
- Airflow Velocity 0.3~0.6m/s
- Pre-filter: Polyester fiber, washable
- LED Display
- Main Body: Cold-rolled steel with anti-bacteria powder coating
- Work Table: 304 stainless steel
- Noise: <60dB
- Front Window: Manual, 5mm toughened glass, anti UV
- Fluorescent Lamp: 8W\*1
- UV Lamp: Emission of 253.7 nanometers
- Consumption: 110W
- Power Supply: 110/220V±10%; 60/50Hz
- Warranty 1 year

### Package Size (WxDxH)

650x550x815mm

### External Size (WxDxH)

550x460x700mm

### Internal Size (WxDxH)

480x340x370mm

**Gross Weight:** 40kg