


FagronLab™ LF-800 UV

Positive pressure workbench



Includes

- 1pc Unguator B/R mixing-homogenizing device.
- 1pc FagronLab™ Pro DMS magnetic Stirrer specially designed for use in Compounding Preparations.
- A set of Compounding Spatulas and Mortar Pestles.



Positive Pressure Airflow Pattern

FagronLab™ LF-800 UV

- Pre-filter: Polyester fiber, washable
- HEPA Filter: 99.999% efficiency at 0.3um
- Front Window: Manual, 5mm toughened glass, anti-UV
- Fluorescent Lamp: 14W*1
- Power Supply: 110/220V±10%; 60/50Hz
- Consumption: 350W
- UV Lamp: Emission of 253.7 nanometers
- Max Opening: 490mm
- Noise: <60dB
- Airflow Velocity: Average of 0.3~0.5m/s, air speed adjustable.
- Main Body: Cold-rolled steel with anti-bacteria powder coating
- Work Table: 304 stainless steel
- Display: LED Display
- Warranty 1 year

Package Size (WxDxH)

960x770x1800mm

External Size (WxDxH)

802x615x1575mm

Internal Size (WxDxH)

800x600x540mm

Work Surface Height: 660mm

Gross Weight: 120kg