


FagronLab™ NP-1 UV

Negative pressure workbench



FagronLab™ NP-1 UV

- HEPA Filter 99.999% efficiency at 0.3um
- Airflow Velocity $\geq 0.5\text{m/s}$
- Pre-filter: Polyester fiber, washable
- LED Display
- Main Body: Cold-rolled steel with anti-bacteria powder coating
- Work Table: 304 stainless steel
- Noise: <60dB
- Front Window: Manual, 5mm toughened glass, anti UV
- Fluorescent Lamp: 8W*1
- UV Lamp: Emission of 253.7 nanometers
- Consumption: 110W
- Power Supply: 110/220V $\pm 10\%$; 60/50Hz
- Warranty 1 year

Package Size (WxDxH)
650x550x815mm

External Size (WxDxH)
550x460x700mm

Internal Size (WxDxH)
480x340x370mm

Gross Weight: 40kg