


FagronLab™ NP-800 UV

Negative pressure workbench



Product Features

Package Size (WxDxH)

106x87x144 cm

External Size (WxDxH)

90x69.5x108 cm

Working area dimensions (WxDxH)

76.8x58x58 cm

Gross Weight

100 kg

- Negative pressure workbench suitable for preparation of volatile odorant products in a secure environment
- Compatible with the capsule filling machine FG1-2
- Special HEPA filtration system with efficiency of 99.999 % at 0.3 um
- Adjustable air circulation volume $\geq 0.5\text{m/s}$
- Function Control via LED display
- Glass sliding front with UV protection filter
- Interior lighting function
- Sterilization function with UV lamp
- Noise <60dB
- Made from stainless steel with antibacterial coating
- Stainless steel work surface
- 1 year warranty

Fagron Hellas

12 km N.R. Trikala - Larisa
P.C. 42100, P.O. Box 32
Trikala, Greece

T +30 24310 83633-5
F +30 24310 83615
www.fagron.gr