

FagronLab™ PRO-DMS

Digital hotplate magnetic stirrer



FagronLab™ PRO-DMS

- Hotplate magnetic stirrer for mixing liquids and preparing solutions
- Quick and easy dissolution of substances with application of temperature
- Electronic control and display of stirring speed and temperature
- Maximum mixing capacity 10.000ml
- Maximum speed 1.500 rpm
- Maximal temperature 340°C
- Warranty 2 years

Package Size (HxWxL)
20x16x34cm

Product Size (HxWxL)
50x16x27cm

Gross Weight:
3kg