

# FagronLab™ RVC1

Laboratory robot vacuum cleaner



## FagronLab™ RVC1

- Professional robot vacuum cleaner for Labs and Health Care Establishments
- Repeatable cleaning mode to guarantee total coverage of laboratory surface and compliance with cleaning SOPs
- UV sterilization function
- Intelligent dust detecting system and mopping function with liquid detergents
- Scheduled cleaning function
- Led display with touch panel
- Extended battery capacity. Vacuuming duration 90 mins
- Auto return to base and charging function
- HEPA filter technology, washable bin
- Remote control for guided cleaning
- Vacuum suction power automatically adjusted to surface and 4 different cleaning modes according to the needs of user (spot, zig-zag, edge and scheduled)
- Loud notifications to inform users about program completion or errors
- Warranty 1 year

### Package Size (LxWxH)

41.2x51.3x15.8cm

### Product Size (DxH)

34x8.5cm

### Gross Weight:

5kg